

INTRODUCING OUR 600 PLATE CAPACITY ANAEROBIC CHAMBER DESIGNED FOR HIGH VOLUME SAMPLES







ANAEROBIC CHAMBER





The BACTRON600 Anaerobic Chamber permits volume handling of anaerobic cultures with maximum operator comfort with 16.0 cu. ft. of workspace up to 300-plate with generous stainless-steel shelving and an additional 300-plate capacity side storage incubator. In combination with workspace capacity and side storage incubator capacity, which is designed for high-volume research requiring nominal bench space than expected for the 600-plate volume plus anaerobic performance features. Our BACTRON Anaerobic Chambers have an airtight construction of stainless steel and rigid acrylic glass for unobstructed vision and integrity. The chamber handcuffs form a comfortable seal around the operator's arms that enables complete manipulation of plates and specimens inside the working chamber forming a work user-friendly in most applications. When sitting undisturbed and unopened at a steady state, the BACTRON Anaerobic Chamber typically rests oxygen-free and anaerobic. Our Anaerobic Chambers meet EN61010, UL61010, and Can/CSA61010 standards.

BACTRONS CREATE AN ENVIRONMENT WITH AN ABSENCE OF OXYGEN WHICH IS OBLIGATE, AND STRICT FOR YOUR ANAEROBIC ORGANISMS

Our BACTRON600 Anaerobic Chamber features:

• Low Gas Supply Alarm Audible and visual notification via the control panel when the AMG gas supply is low. Also, operators can change the cylinder while the chamber atmosphere stays anaerobic.

• Auto-Sleeve Cycle No more guesswork or operator variation; automated ultra-efficient sleeve cycles with one press of the foot switch.

Integrated Internal Regulators

Ensure consistent and reproducible gas delivery to all chamber processes plus protecting all the internal components from incorrect cylinder settings.

User Adjustable Airlock
Cycling

Customize your airlock cycle to fit your time and budget allowing operator customization from 3, 4, or 5 auto-cycle iterations.

• Auto-Commissioning Cycle

The single-button autocommissioning cycle automatically creates an anaerobic atmosphere in 5 hours; press the button and walk away while the chamber becomes anaerobic.

FOLLOW US ON:





HELDO

PN: 0740565 1/23

0-	INCUBATOR	
STATUS	SET TEMPERATURE	
	WORKSPACE	
AIRLOCK		
	SLEEVES	
AUTO CYCLE AIRLOCK PRESSURE	SLEEVE CYCLE IN PROGRESS ((()) RIGH FOOT REPLAY TO START RIGH ACAR TO ABORT CONTRIBUTION OF THE ADORT	
	SLEEVE CYCLE SLEEVE PRESSURE	

Applications:

- » Clinical Microbiology
- » Microbiome Studies
- » Susceptibility Assays
- » Fermentation and Spoilage Studies
- » Food Microbiology
- » Bio-Energy Research
- » Environmental Microbiology
- » Anaerobic Bacteriology

Accessories		Unit Specifications	
Item Description	Part Number	BACTRON600	Part ID: BAA60022/ [BAA60022-E]*
Handheld Oxygen Sensor	9902223	Unit Weight	412lbs/187kg
Activated Carbon Media (2lbs/0.9kgs)	1060500	Shipping Weight	657lbs/298kg
Petri Plate Racks (2 x 11)	5110729	Workspace Chamber Volume	16.0 cu.ft. /453L
Leak Detector	4600501	Workspace Chamber Dimensions	42.5" x 28.9" x 25.0" 1059 x 734 x 635mm
UV Viewing Lamp	9490507	- W x D x H	
Nitrogen Regulator Kit	9740567	Exterior Dimensions W x D x H	88.5" x 32.5" x 27.6" 2248 x 826 x 701mm
Nitrile Cuffs	3600525(S) 3600526(M) 3600527(L)	Airlock Dimensions W x D x Diagonal	16.0" x 10.0" x 11.5" 406 x 254 x 292mm
Activated Carbon, Volatile Com-	9490578 (110V)	Airlock Volume	1.29 cu.ft. /36.5L
pounds Scrubber Fan	9490581 (220V)	Airlock Plate Capacity	216
BACTRON600 Rolling Stand		Incubator Specifications	
(W x D x H) 88.5" x 30" x 29"/2250 x 762 x 762 mm	BACSTAND-LG22	Plate Capacity	600
AMG REGULATOR ASSY 3/16ID		Workspace Incubator Volume	4.6 cu.ft. /130L
TUBING / A HIGHER-PRESSURE GAUGE W/WIDE TUBING	9740560	Side Incubator Dimensions	23.5" x 18.5" 597 x 470mm
Stereo Microscope w/ Accessories	Contact <u>sales@sheldonmfg.</u> <u>com</u>	Temperature Range	Ambient + 5° to 70° C
Anaerobic Chamber Start-Up Kit*	9490511	Temperature Uniformity	+/- 1.0° C @ 37°C
	3600500(S)	Side Incubator Uniformity	+/- 1.0° C @ 37°C
Latex Cuffs	3600502(L)	Power	
Benzalkonium Chloride Solution	1060501	AC Voltage	110-120/220-240
Chamber Cleaner		Amperage	11/8
Anaerobic Indicator Strips	1060553	Frequency	50-60Hz
Activated Charcoal Media	1060500	*Denotes 220v unit	
Airlock Gasket (12 x 12)	3450507		
Acrylic Glass Cleaner	1060503		
Acrylic Glass Scratch Remover	1060504		
O-Ring For Sleeve Cuff	6000504		
	1		



https://www.sheldonmanufacturing.com

sales@sheldonmfg.com

Manufactured in Cornelius, Oregon USA



SHELD

USER MANUAI

0

PN: 0740565 1/23

FOLLOW US ON: Pin f



DESIGNED FOR HIGHER VOLUME RESEARCH

