

Precise Temperature Control - Superior Uniformity

- Independent Over Temperature Thermostat
- Over Temperature Alarm
- Temperature Range 0°C to 45°C at 20°C Ambient



SRI3



		Model Number
Refrigerated Incubators	110-120V	SRI3
	220-230V	SRI3-2
Details		Under Counter
Exterior Dimensions (wxdxh)	Inches	24 x 21 x 33.8
	N m	610 x 534 x 858
Chamber Dimensions (wxdxh)	Inches	16 x 1 x 2.5
	N m	40 x 3 x 5
Incubator Chamber Capacity	cu ft	2.
	L	6
Interior Outlet	Number	110V / 220V
Temperature Uniformity	Celsius	+/- 0.5°C at 20°C
Bottle Capacity	Number	62
Number of Shelves	Included	2

*All specifications are determined by using average values on standard equipment at an ambient temperature of 25°C (77°F) and line voltages within +/-10% of unit type (115V/230V). Temperature specifications follow DIN 12880 methodology. We reserve the right to change specifications at any time.
 For a complete list of product specifications, call (800)322-4897 or visit the SHEL LAB website, www.TIFMEPONGH.com

0740506 /