

Lab ARMOR® Tall Vessel Lid

Introducing Lab Armor Tall Vessel Lids! With the added height, these new Lab Armor Tall Vessel Lids offer room for taller containers and provide more protection for samples.

An easily attachable design, the Lab Armor Tall Vessel Lids install in any corresponding Lab Armor Bead Bath.



Cell culture in bottles requires special attention to prevent contamination of water. Lab Armor Tall Vessel Lids are the perfect height to accommodate cell culture in culture flasks.

Features:

- Stainless steel
- Hinged style
- Easy installation
- Taller lid provides room for large cell culture flasks
- Available in three sizes for 6 L, 14 L, 20 L baths



	6L Bead Bath Lid	14L Bead Bath Lid	20L Bead Bath Lid
Part ID/SKU	9521387	9521388	9521389
Dimensions (w x d x h)	12.6" x 7.5" x 5" (320 x 191 x 127 mm)	12.6" x 13.5" x 5" (320 x 343 x 127 mm)	12.6" x 18.5" x 5" (320 x 470 x 127mm)
Overhead Clearance	16.5" (419 mm)	22.5" (572 mm)	27.5" (699 mm)

Order a Lab Armor Tall Vessel Lid Today!

sales@labarmor.com | 1-888-227-1410 | labarmor.com